

# Crono P Pulsatile Infusion Pump

**For subcutaneous infusion applications.**

The Crono P is a new pump designed to provide a pulsatile drug delivery pattern. It incorporates a real time clock and has the ability to set three different pulse rates for cyclical delivery patterns. It has a basic setting which allows a choice of 90 or 180 minute intervals to suit a wide range of applications.

The pump uses dedicated 10 ml or 20 ml syringe reservoirs, and is supplied with a range of carrying accessories for user convenience. The single battery cell has a typical life expectancy sufficient for 200-300 full syringe infusions.



## Specifications

Dimensions	76 x 49 x 29 mm
Weight	115 g (Battery Included)
Battery	3 Volt Type 123A Lithium
Syringe Volume	Select 10 or 20ml
Volume Administration	Pulsatile pattern with three programmable levels using real time clock for cyclical infusions.
Accuracy	+/- 2%
Keyboard Lock	2 Security levels available

## Accessories

Each pump is supplied with a range of accessories for patient convenience, which include:

- Soft Case
- Battery cover removal tool
- Soft waist belt
- Neck Strap
- Spare battery

LR004/19BM