


# Quick-Set™ Infusion Set

The Quick-Set™ infusion set is suitable for all long term applications of subcutaneous drug infusion, and is widely used for the administration of insulin, diamorphine and a range of other drugs. The Quick-Set™ has a vertically mounted catheter to ensure consistent depth of penetration, and offers two catheter lengths to suit a wide variety of patient types. It also incorporates an at-site detachable connector. An insertion aid, the Qickserter, is available with these devices (Product code MMT 305QS)

Code	Description	NHS Catalogue	Catheter Length	Tubing Length
MMT-390	Quick Set	FSB1899	9mm	110cm
MMT-391	Quick Set	FSB1900	6mm	110cm
MMT-392	Quick Set	FSB1901	9mm	60cm
MMT-393	Quick Set	FSB1902	6mm	60cm


Straight-in insertion with 6mm or 9mm options to suit a variety of patient types and applications



Soft polyurethane catheter with disconnectable tubing set for comfort and convenience



Alignment indicator to unlock or remove connector



Very low volume kink-free tubing holds less than 0.2ml. No PVC contact with fluid path for drug compatibility

Integral fixing tape for convenience

Low profile connector system

Optional inserter device for patient



LR008/19BM