

neria™ Infusion Set

For the highest levels of comfort, security and safety in subcutaneous infusion practice

The neria™ infusion set is suitable for all applications of subcutaneous drug infusion, and is widely used for the administration of diamorphine and a range of other drugs. The fixed needle length and orientation ensure controlled depth of penetration. The set is available in the following configurations. (Formerly known as Contact or Thalaset)

AVAILABLE ON PRESCRIPTION

Code	Description	Needle Gauge	Needle Length	Tubing Length	PIP Code	NPC Code
78-060-2738	Neria G27	27g	8mm	60cm	N/A	FSN133
78-110-2738	Neria G27	27g	8mm	110cm	353-0979	FSN127
78-110-2731	Neria G27	27g	10mm	110cm	353-0946	FSN128
78-060-2938	Neria G29	29g	8mm	60cm	353-0987	FSN129
78-060-2931	Neria G29	29g	10mm	60cm	353-0995	FSN131
78-110-2938	Neria G29	29g	8mm	110cm	353-1019	FSN130
78-110-2931	Neria G29	29g	10mm	110cm	353-1035	FSN132



LR009/19BM