


neria™ multi Infusion Set

For simplified delivery of higher volumes in immunoglobulin infusion therapies

Based on the well established neria™ line of subcutaneous infusion sets, the new neria multi range addresses the requirements for short term, high volume infusions, primarily for immunoglobulin administration.

With 27g needles, microvolume tubing sets and a choice of 2, 3 or 4 way configurations these sets should address the needs of the widest range of patients. The low volume characteristics will offer simple priming procedures and significant reduction in drug wastage.



Code	Description	Needle Gauge	Needle Length	Tubing Length
721060-2738	2 way set	27g	8 mm	100cm
721060-2731	2 way set	27g	10 mm	100cm
721060-3738	3 way set	27g	8 mm	100cm
721060-3731	3 way set	27g	10 mm	100cm
721060-3737	3 way set	27g	12 mm	100cm
721060-4738	4 way set	27g	8 mm	100cm
721060-4731	4 way set	27g	10 mm	100cm
721060-4737	4 way set	27g	12 mm	100cm



LR011/19BM