

# CRONOJET 100



For fast and easy filling  
of the 100 ml crn<sup>®</sup> CRONO<sup>®</sup> Reservoir

# CRONOJET 100



**SCALE 1:1**

## DESCRIPTION

Device designed for fast and easy filling of the 100 ml crn<sup>®</sup> CRONO<sup>®</sup> Reservoir using a spike with venting filter.

## CLASSIFICATION

Class I medical device 

• **SUPPLY VOLTAGE**  
3V dc - 2 AA batteries

• **DIMENSIONS**  
124 x 71 x 56 mm

• **WEIGHT**  
178 g batteries included

**Manufactured by:**



**CANÈ S.p.A.**  
via Cuorgnè, 42/a - 10098 Rivoli (TO) - Italy  
Tel. 011-9574872 - Fax 011-9598880  
[www.canespa.it](http://www.canespa.it) - [mailbox@canespa.it](mailto:mailbox@canespa.it)

**Distributed by:**



**Applied Medical Technology Ltd**  
3960 Enterprise, Cambridge Research Park, Waterbeach CB25 9PE  
Tel: 01223 420415 • e-mail: [sales@applied-medical.co.uk](mailto:sales@applied-medical.co.uk)  
[www.applied-medical.co.uk](http://www.applied-medical.co.uk)